

Product datasheet for **TP760649**

PDE6D (NM_002601) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human phosphodiesterase 6D, cGMP-specific, rod, delta (PDE6D), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length PDE6D
Tag:	N-His
Predicted MW:	17.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	20 mM PB, pH 7.4
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_002592
Locus ID:	5147
UniProt ID:	O43924 , Q6IB24
RefSeq Size:	1214
Cytogenetics:	2q37.1
RefSeq ORF:	450
Synonyms:	JBTS22; PDED



[View online »](#)

Summary:

This gene encodes the delta subunit of rod-specific photoreceptor phosphodiesterase (PDE), a key enzyme in the phototransduction cascade. A similar protein in cow functions in solubilizing membrane-bound PDE. In addition to its role in the PDE complex, the encoded protein is thought to bind to prenyl groups of proteins to target them to subcellular organelles called cilia. Mutations in this gene are associated with Joubert syndrome-22. Alternative splicing results in multiple splice variants. [provided by RefSeq, Mar 2014]

Protein Pathways:

Progesterone-mediated oocyte maturation, Purine metabolism

Product images: