

Product datasheet for **TP760711**

CBR4 (NM_032783) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human carbonyl reductase 4 (CBR4), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length CBR4
Tag:	N-His
Predicted MW:	25.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_116172
Locus ID:	84869
UniProt ID:	Q8N4T8
RefSeq Size:	3523
Cytogenetics:	4q32.3
RefSeq ORF:	711
Synonyms:	SDR45C1



[View online »](#)

Summary:

The heterotetramer with HSD17B8 has NADH-dependent 3-ketoacyl-acyl carrier protein reductase activity, and thereby plays a role in mitochondrial fatty acid biosynthesis (PubMed:19571038, PubMed:25203508). Within the heterotetramer, HSD17B8 binds NADH; CBR4 binds NADPD (PubMed:25203508). The homotetramer has NADPH-dependent quinone reductase activity (PubMed:19000905). Both homotetramer and the heterotetramer have broad substrate specificity and can reduce 9,10-phenanthrenequinone, 1,4-benzoquinone and various other o-quinones and p-quinones (in vitro) (PubMed:19000905, PubMed:19571038, PubMed:25203508).[UniProtKB/Swiss-Prot Function]

Protein Families:

Druggable Genome

Product images: