

Product datasheet for **TP760724**

FAM54A (MTRF2) (NM_001099286) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human family with sequence similarity 54, member A (FAM54A), transcript variant 1, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length FAM54A
Tag:	N-His
Predicted MW:	43.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001092756
Locus ID:	113115
UniProt ID:	Q6P444
RefSeq Size:	1811
Cytogenetics:	6q23.3
RefSeq ORF:	1155
Synonyms:	DUFD1; FAM54A
Summary:	May play a role in mitochondrial aerobic respiration essentially in the testis. Can also promote mitochondrial fission (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

