

Product datasheet for **TP760742**

ETV3 (NM_005240) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human ets variant 3 (ETV3), transcript variant 2, full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length ETV3
Tag:	N-His
Predicted MW:	16.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005231
Locus ID:	2117
UniProt ID:	P41162
RefSeq Size:	1443
Cytogenetics:	1q23.1
RefSeq ORF:	429
Synonyms:	METS; PE-1; PE1
Summary:	Transcriptional repressor that contribute to growth arrest during terminal macrophage differentiation by repressing target genes involved in Ras-dependent proliferation. Represses MMP1 promoter activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Protein Families: Transcription Factors

Product images:

