

Product datasheet for **TP760864**

TBLR1 (TBL1XR1) (NM_024665) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human transducin (beta)-like 1 X-linked receptor 1 (TBL1XR1), full length, with N-terminal HIS tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Tag:	N-His
Predicted MW:	55.4 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25mM Tris, pH8.0, 150mM NaCl, 10% glycerol, 1 % Sarkosyl.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_078941
Locus ID:	79718
RefSeq Size:	6550
Cytogenetics:	3q26.32
RefSeq ORF:	1542
Synonyms:	C21; DC42; IRA1; MRD41; TBLR1
Summary:	This gene is a member of the WD40 repeat-containing gene family and shares sequence similarity with transducin (beta)-like 1X-linked (TBL1X). The protein encoded by this gene is thought to be a component of both nuclear receptor corepressor (N-CoR) and histone deacetylase 3 (HDAC 3) complexes, and is required for transcriptional activation by a variety of transcription factors. Mutations in these gene have been associated with some autism spectrum disorders, and one finding suggests that haploinsufficiency of this gene may be a cause of intellectual disability with dysmorphism. Mutations in this gene as well as recurrent translocations involving this gene have also been observed in some tumors. [provided by RefSeq, Mar 2016]



[View online »](#)

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Wnt signaling pathway

Product images:

