

Product datasheet for TP760893

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

RNASE9 (NM_001110360) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human ribonuclease, RNase A family, 9 (non-active)

(RNASE9), transcript variant 3, full length, with N-terminal HIS tag, expressed in E. coli, 50ug

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

A DNA sequence encoding human full-length RNASE9

Tag: N-His

Predicted MW: 24.6 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 50 mM Tris-HCl, pH 8.0, 500 mM NaCl, 10% glycerol

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeg: NP 001103830

 Locus ID:
 390443

 UniProt ID:
 P60153

 Cytogenetics:
 14q11.2

 RefSeq ORF:
 630

Synonyms: h461; HEL128; RAK1

Summary: Does not exhibit any ribonuclease activity.[UniProtKB/Swiss-Prot Function]

Protein Families: Secreted Protein





Product images:

