

Product datasheet for **TP760959**

ZNF652 (NM_014897) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human zinc finger protein 652 (ZNF652), transcript variant 2, full length, with N-terminal HIS tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length ZNF652
Tag:	N-His
Predicted MW:	69.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 8 M urea
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_055712
Locus ID:	22834
UniProt ID:	Q9Y2D9
RefSeq Size:	5428
Cytogenetics:	17q21.32-q21.33
RefSeq ORF:	1818
Summary:	Functions as a transcriptional repressor.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

