

Product datasheet for **TP760974**

C2orf84 (FAM228A) (NM_001040710) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human chromosome 2 open reading frame 84 (C2orf84), full length, with N-terminal HIS tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length C2orf84
Tag:	N-His
Predicted MW:	23.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 8 M urea
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001035800</u>
Locus ID:	653140
UniProt ID:	<u>Q86W67</u>
RefSeq Size:	1889
Cytogenetics:	2p23.3
RefSeq ORF:	618
Synonyms:	C2orf84


[View online »](#)

Product images:

