

## Product datasheet for **TP760995**

### MPZL3 (NM\_198275) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human myelin protein zero-like 3 (MPZL3), full length, with N-terminal HIS tag, expressed in E.coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length MPZL3
Tag:	N-His
Predicted MW:	25.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 8 M urea
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_938016</a>
Locus ID:	196264
UniProt ID:	<a href="#">Q6UWV2</a>
RefSeq Size:	3986
Cytogenetics:	11q23.3
RefSeq ORF:	705
Summary:	Mediates homophilic cell-cell adhesion.[UniProtKB/Swiss-Prot Function]
Protein Families:	Transmembrane



[View online »](#)

**Product images:**

