

Product datasheet for **TP761050**

BTBD19 (NM_001136537) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human BTB (POZ) domain containing 19 (BTBD19), full length, with N-terminal HIS tag, expressed in E.coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length BTBD19
Tag:	N-His
Predicted MW:	32.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001130009</u>
Locus ID:	149478
UniProt ID:	<u>C9JJ37</u>
Cytogenetics:	1p34.1
RefSeq ORF:	873

[View online »](#)

Product images:

