

Product datasheet for **TP761063**

PDZD4 (NM_032512) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human PDZ domain containing 4 (PDZD4), full length, with N-terminal HIS tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length PDZD4
Tag:	N-His
Predicted MW:	86 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_115901
Locus ID:	57595
UniProt ID:	Q76G19
RefSeq Size:	3745
Cytogenetics:	Xq28
RefSeq ORF:	2307
Synonyms:	LNx5; LU1; PDZK4; PDZRN4L



[View online »](#)

Product images:

