

## Product datasheet for **TP761075**

### **FLJ14213 (PRR5L) (NM\_024841) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Human proline rich 5 like (PRR5L), transcript variant 2, full length, with N-terminal HIS tag, expressed in E. coli, 50ug
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	A DNA sequence encoding human full-length PRR5L
<b>Tag:</b>	N-His
<b>Predicted MW:</b>	40.7 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_079117</a>
<b>Locus ID:</b>	79899
<b>UniProt ID:</b>	<a href="#">Q6MZQ0</a>
<b>RefSeq Size:</b>	3946
<b>Cytogenetics:</b>	11p13-p12
<b>RefSeq ORF:</b>	1104
<b>Synonyms:</b>	PROTOR2



[View online »](#)

**Summary:**

Associates with the mTORC2 complex that regulates cellular processes including survival and organization of the cytoskeleton (PubMed:17461779). Regulates the activity of the mTORC2 complex in a substrate-specific manner preventing for instance the specific phosphorylation of PKCs and thereby controlling cell migration (PubMed:22609986). Plays a role in the stimulation of ZFP36-mediated mRNA decay of several ZFP36-associated mRNAs, such as TNF-alpha and GM-CSF, in response to stress (PubMed:21964062). Required for ZFP36 localization to cytoplasmic stress granule (SG) and P-body (PB) in response to stress (PubMed:21964062).[UniProtKB/Swiss-Prot Function]

**Product images:**