

## Product datasheet for **TP761116**

### **C10orf84 (FAM204A) (NM\_022063) Human Recombinant Protein**

#### **Product data:**

|  |  |
|--|--|
| <b>Product Type:</b>                         | Recombinant Proteins   |
| <b>Description:</b>                          | Purified recombinant protein of Human chromosome 10 open reading frame 84 (C10orf84), transcript variant 1, full length, with N-terminal HIS tag, expressed in E. coli, 50ug |
| <b>Species:</b>                              | Human  |
| <b>Expression Host:</b>                      | E. coli  |
| <b>Expression cDNA Clone or AA Sequence:</b> | A DNA sequence encoding human full-length C10orf84   |
| <b>Tag:</b>                                  | N-His  |
| <b>Predicted MW:</b>                         | 26.8 kDa   |
| <b>Concentration:</b>                        | >0.05 µg/µL as determined by microplate BCA method   |
| <b>Purity:</b>                               | > 80% as determined by SDS-PAGE and Coomassie blue staining  |
| <b>Buffer:</b>                               | 25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol   |
| <b>Note:</b>                                 | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.                         |
| <b>Storage:</b>                              | Store at -80°C.  |
| <b>Stability:</b>                            | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.                                |
| <b>RefSeq:</b>                               | <a href="#">NP_071346</a>  |
| <b>Locus ID:</b>                             | 63877  |
| <b>UniProt ID:</b>                           | <a href="#">Q9H8W3</a>   |
| <b>RefSeq Size:</b>                          | 2766   |
| <b>Cytogenetics:</b>                         | 10q26.11   |
| <b>RefSeq ORF:</b>                           | 699  |
| <b>Synonyms:</b>                             | bA319I23.1; C10orf84   |



[View online »](#)

**Product images:**

