

## Product datasheet for **TP761202**

### **ACTR6 (NM\_022496) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Human ARP6 actin-related protein 6 homolog (yeast) (ACTR6), full length, with N-terminal HIS tag, expressed in E. coli, 50ug
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	A DNA sequence encoding human full-length ACTR6
<b>Tag:</b>	N-His
<b>Predicted MW:</b>	45.6 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	50 mM Tris-HCl, pH 8.0, 8 M urea
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_071941</a>
<b>Locus ID:</b>	64431
<b>UniProt ID:</b>	<a href="#">Q9GZN1</a>
<b>RefSeq Size:</b>	2491
<b>Cytogenetics:</b>	12q23.1
<b>RefSeq ORF:</b>	1188
<b>Synonyms:</b>	ARP6; CDA12; hARP6; hARPX; HSPC281; MSTP136



[View online »](#)

**Product images:**

