

Product datasheet for TP761273

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ZNF133 (NM_001083330) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human zinc finger protein 133 (ZNF133), transcript variant 2,

full length, with N-terminal HIS tag, expressed in E. coli, 50ug

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

A DNA sequence encoding human full-length ZNF133

Tag: N-His

Predicted MW: 73.2 kDa

Concentration: $>0.05 \mu g/\mu L$ as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 50 mM Tris-HCl, pH 8.0, 8 M urea

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeg: NP 001076799

 Locus ID:
 7692

 UniProt ID:
 P52736

 RefSeq Size:
 2500

 Cytogenetics:
 20p11.23

RefSeq ORF: 1959

Synonyms: pHZ-13; pHZ-66; ZNF150

Summary: May be involved in transcriptional regulation as a repressor.[UniProtKB/Swiss-Prot Function]

Protein Families: Transcription Factors





Product images:

