

## Product datasheet for **TP761470**

### **PPIAL4E (NM\_001144032) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Human peptidylprolyl isomerase A (cyclophilin A)-like 4E (PPIAL4E), full length, with N-terminal HIS tag, expressed in E.coli, 50ug
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	A DNA sequence encoding human full-length PPIAL4E
<b>Tag:</b>	N-His
<b>Predicted MW:</b>	18 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_001137504</a>
<b>Locus ID:</b>	730262
<b>UniProt ID:</b>	<a href="#">A0A075B759</a> , <a href="#">P0DN26</a>
<b>Cytogenetics:</b>	1q21.1
<b>RefSeq ORF:</b>	492
<b>Synonyms:</b>	COAS2; PPIAL4F
<b>Summary:</b>	PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Product images:**

