

Product datasheet for **TP761554**

MAB21L3 (NM_152367) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human mab-21-like 3 (<i>C. elegans</i>) (MAB21L3), full length, with N-terminal HIS tag, expressed in <i>E. coli</i> , 50ug
Species:	Human
Expression Host:	<i>E. coli</i>
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length MAB21L3
Tag:	N-His
Predicted MW:	42.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_689580
Locus ID:	126868
UniProt ID:	Q8N8X9
RefSeq Size:	3229
Cytogenetics:	1p13.1
RefSeq ORF:	1086
Synonyms:	C1orf161



[View online »](#)

Product images: