

Product datasheet for **TP761601**

GASK1B (NM_001031700) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human family with sequence similarity 198, member B (FAM198B), transcript variant 1, full length, with N-terminal GST and C-terminal HIS tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length FAM198B
Tag:	N-GST and C-His
Predicted MW:	86.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001026870
Locus ID:	51313
UniProt ID:	Q6UWH4 , Q9NYZ0
RefSeq Size:	4833
Cytogenetics:	4q32.1
RefSeq ORF:	1581
Synonyms:	AD021; AD036; C4orf18; ENED; FAM198B



[View online »](#)

Product images:

