

Product datasheet for **TP761625**

RCOR3 (NM_018254) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human REST corepressor 3 (RCOR3), transcript variant 4, full length, with N-terminal GST and C-terminal HIS tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length RCOR3
Tag:	N-GST and C-His
Predicted MW:	83.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_060724
Locus ID:	55758
UniProt ID:	Q9P2K3
RefSeq Size:	4083
Cytogenetics:	1q32.2-q32.3
RefSeq ORF:	1485
Summary:	May act as a component of a corepressor complex that represses transcription. [UniProtKB/Swiss-Prot Function]
Protein Families:	Transcription Factors



[View online »](#)

Product images:

