

## Product datasheet for **TP761673**

### **SAMD10 (NM\_080621) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Human sterile alpha motif domain containing 10 (SAMD10), full length, with N-terminal HIS tag, expressed in E. coli, 50ug
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	A DNA sequence encoding human full-length SAMD10
<b>Tag:</b>	N-His
<b>Predicted MW:</b>	22.6 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	50 mM Tris-HCl, pH 8.0, 8 M urea
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_542188</a>
<b>Locus ID:</b>	140700
<b>UniProt ID:</b>	<a href="#">Q9BYL1</a>
<b>RefSeq Size:</b>	2196
<b>Cytogenetics:</b>	20q13.33
<b>RefSeq ORF:</b>	606
<b>Synonyms:</b>	C20orf136; dj591C20; dj591C20.7



[View online »](#)

**Product images:**