

Product datasheet for **TP761681**

TOMM40 (NM_001128916) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human translocase of outer mitochondrial membrane 40 homolog (yeast) (TOMM40), nuclear gene encoding mitochondrial protein, transcript variant 2, full length, with N-terminal HIS tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Tag:	N-His
Predicted MW:	37.7 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50mM Tris, pH8.0,8M Urea.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001122388
Locus ID:	10452
UniProt ID:	O96008
Cytogenetics:	19q13.32
RefSeq ORF:	1083
Synonyms:	C19orf1; D19S1177E; PER-EC1; PEREC1; TOM40
Summary:	The protein encoded by this gene is localized in the outer membrane of the mitochondria. It is the channel-forming subunit of the translocase of the mitochondrial outer membrane (TOM) complex that is essential for import of protein precursors into mitochondria. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Aug 2015]
Protein Families:	Druggable Genome, Ion Channels: Other
Protein Pathways:	Amyotrophic lateral sclerosis (ALS)



[View online »](#)

Product images: