

Product datasheet for **TP761695**

LBX2 (NM_001009812) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human ladybird homeobox 2 (LBX2), full length, with N-terminal HIS tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length LBX2
Tag:	N-His
Predicted MW:	21.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 8 M urea
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001009812
Locus ID:	85474
UniProt ID:	Q6XYB7
RefSeq Size:	1463
Cytogenetics:	2p13.1
RefSeq ORF:	582
Synonyms:	LP3727
Summary:	Putative transcription factor.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

