

## Product datasheet for **TP761778**

### **DULLARD (CTDNEP1) (NM\_015343) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Human CTD nuclear envelope phosphatase 1 (CTDNEP1), transcript variant 1, full length, with N-terminal GST and C-terminal His tag, expressed in E. coli, 50ug
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	A DNA sequence encoding human full-length CTDNEP1
<b>Tag:</b>	N-GST and C-His
<b>Predicted MW:</b>	56.2 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	50 mM Tris-HCl, pH 8.0, 8 M urea
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_056158</a>
<b>Locus ID:</b>	23399
<b>UniProt ID:</b>	<a href="#">O95476</a>
<b>RefSeq Size:</b>	1780
<b>Cytogenetics:</b>	17p13.1
<b>RefSeq ORF:</b>	732
<b>Synonyms:</b>	DULLARD; HSA011916; NET56



[View online »](#)

**Summary:**

Serine/threonine protein phosphatase forming with CNEP1R1 an active phosphatase complex that dephosphorylates and may activate LPIN1 and LPIN2. LPIN1 and LPIN2 are phosphatidate phosphatases that catalyze the conversion of phosphatidic acid to diacylglycerol and control the metabolism of fatty acids at different levels. May indirectly modulate the lipid composition of nuclear and/or endoplasmic reticulum membranes and be required for proper nuclear membrane morphology and/or dynamics. May also indirectly regulate the production of lipid droplets and triacylglycerol. May antagonize BMP signaling. [UniProtKB/Swiss-Prot Function]

**Protein Families:**

Transmembrane

**Product images:**