

Product datasheet for **TP761860**

PPIAL4A (NM_178230) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human peptidylprolyl isomerase A (cyclophilin A)-like 4A (PPIAL4A), full length, with N-terminal GST and C-terminal His tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length PPIAL4A
Tag:	N-GST and C-His
Predicted MW:	46 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 8 M urea
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_839944
Locus ID:	164022
RefSeq Size:	660
Cytogenetics:	1q21.2
RefSeq ORF:	492
Synonyms:	COAS-2; COAS2; PPIAL4


[View online »](#)

Product images:

