

Product datasheet for **TP761894**

C1orf54 (NM_024579) Human Recombinant Protein

Product data:

| | |
|--|---|
| Product Type: | Recombinant Proteins |
| Description: | Purified recombinant protein of Human chromosome 1 open reading frame 54 (C1orf54), full length, with N-terminal GST and C-terminal His tag, expressed in E. coli, 50ug |
| Species: | Human |
| Expression Host: | E. coli |
| Expression cDNA Clone or AA Sequence: | A DNA sequence encoding human full-length C1ORF54 |
| Tag: | N-GST and C-His |
| Predicted MW: | 42.7 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_078855 |
| Locus ID: | 79630 |
| UniProt ID: | Q8WWF1 |
| RefSeq Size: | 543 |
| Cytogenetics: | 1q21.2 |
| RefSeq ORF: | 393 |
| Protein Families: | Secreted Protein |



[View online »](#)

Product images:

