

## Product datasheet for **TP761956**

### HSZFP36 (ZNF823) (NM\_001080493) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of ZNF823,full length, with N-terminal GST and C-terminal His tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding full-length of ZNF823
Tag:	N-GST and C-His
Predicted MW:	98.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001073962</a>
Locus ID:	55552
UniProt ID:	<a href="#">P16415</a>
RefSeq Size:	2478
Cytogenetics:	19p13.2
RefSeq ORF:	1830
Synonyms:	HSZFP36
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Product images:**