

Product datasheet for **TP761986**

SLC30A4 (NM_013309) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Purified recombinant protein of Human solute carrier family 30 (zinc transporter), member 4 (SLC30A4),Pro336-Gln425, with N-terminal His-ABP tag, expressed in E. coli, 50ug |
| Species: | Human |
| Expression Host: | E. coli |
| Expression cDNA Clone or AA Sequence: | A DNA sequence encoding the region(Pro336-Gln425)of SLC30A4 |
| Tag: | N-His-ABP (Albumin-Binding Protein) |
| Predicted MW: | 25.3 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 10% glycerol |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_037441 |
| Locus ID: | 7782 |
| UniProt ID: | O14863 |
| RefSeq Size: | 4290 |
| Cytogenetics: | 15q21.1 |
| RefSeq ORF: | 1287 |
| Synonyms: | znT-4; ZNT4 |



[View online »](#)

Summary:

Zinc is the second most abundant trace metal in the human body. It is an essential element, serving both a structural role, as in the formation of zinc fingers in DNA-binding proteins, and a catalytic role in metalloenzymes, such as pancreatic carboxypeptidases (e.g., MIM 114852), alkaline phosphatases (e.g., MIM 171760), various dehydrogenases, and superoxide dismutases (e.g., MIM 147450). SLC30A4, or ZNT4, belongs to the ZNT family of zinc transporters. ZNTs are involved in transporting zinc out of the cytoplasm and have similar structures, consisting of 6 transmembrane domains and a histidine-rich cytoplasmic loop (Huang and Gitschier, 1997 [PubMed 9354792]).[supplied by OMIM, Mar 2008]

Protein Families:

Transmembrane

Product images: