

Product datasheet for **TP762085**

CDw75 (ST6GAL1) (NM_173216) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human ST6 beta-galactosamide alpha-2,6-sialyltransferase 1 (ST6GAL1), transcript variant 1,Lys27-End, with N-terminal His tag, expressed in E. coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding the region(Lys27-End) of ST6GAL1
Tag:	N-His
Predicted MW:	43.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 8 M urea
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_775323
Locus ID:	6480
UniProt ID:	P15907 , B2R513
RefSeq Size:	4604
Cytogenetics:	3q27.3
RefSeq ORF:	1218
Synonyms:	SIAT1; ST6Gall; ST6N



[View online »](#)

Summary:

This gene encodes a member of glycosyltransferase family 29. The encoded protein is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The protein, which is normally found in the Golgi but can be proteolytically processed to a soluble form, is involved in the generation of the cell-surface carbohydrate determinants and differentiation antigens HB-6, CD75, and CD76. This gene has been incorrectly referred to as CD75. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2017]

Protein Families:

Secreted Protein

Protein Pathways:

Metabolic pathways, N-Glycan biosynthesis

Product images: