

Product datasheet for **TP762190**

ANGPTL1 (NM_004673) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human angiotensin-like 1 (ANGPTL1), Gly24-End, with N-terminal His tag, expressed in E.coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding the region(Gly24-End) of ANGPTL1
Tag:	N-His
Predicted MW:	54.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 8 M urea
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004664
Locus ID:	9068
UniProt ID:	O95841
RefSeq Size:	3555
Cytogenetics:	1q25.2
RefSeq ORF:	1473
Synonyms:	ANG3; ANGPT3; AngY; ARP1; dj595C2.2; UNQ162




[View online »](#)

Summary:

Angiopoietins are members of the vascular endothelial growth factor family and the only known growth factors largely specific for vascular endothelium. Angiopoietin-1, angiopoietin-2, and angiopoietin-4 participate in the formation of blood vessels. The protein encoded by this gene is another member of the angiopoietin family that is widely expressed in adult tissues with mRNA levels highest in highly vascularized tissues. This protein was found to be a secretory protein that does not act as an endothelial cell mitogen in vitro. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome, Secreted Protein

Product images:

Purified recombinant protein ANGPTL1 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.