

Product datasheet for **TP762363**

TdT (DNNT) (NM_001017520) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human deoxynucleotidyltransferase, terminal (DNNT), transcript variant 2, Met1-Thr160, with N-terminal His tag, expressed in E.coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding the region(Met1-Thr160) of TDT
Tag:	N-His
Predicted MW:	18.0 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001017520
Locus ID:	1791
UniProt ID:	P04053
RefSeq Size:	2068
Cytogenetics:	10q24.1
RefSeq ORF:	1527
Synonyms:	TDT




[View online »](#)

Summary:

This gene is a member of the DNA polymerase type-X family and encodes a template-independent DNA polymerase that catalyzes the addition of deoxynucleotides to the 3'-hydroxyl terminus of oligonucleotide primers. In vivo, the encoded protein is expressed in a restricted population of normal and malignant pre-B and pre-T lymphocytes during early differentiation, where it generates antigen receptor diversity by synthesizing non-germ line elements (N-regions) at the junctions of rearranged Ig heavy chain and T cell receptor gene segments. Alternatively spliced transcript variants encoding different isoforms of this gene have been described. [provided by RefSeq, Jul 2008]

Protein Pathways:

Hematopoietic cell lineage, Non-homologous end-joining

Product images:

Purified recombinant protein TDT was analyzed by SDS-PAGE gel and Coomassie Blue Staining.