

Product datasheet for **TP762452**

NSE (ENO2) (NM_001975) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human enolase 2 (gamma, neuronal) (ENO2), Val188-Glu293, with N-terminal His tag, expressed in E.coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding the region(Val188-Glu293) of ENO2
Tag:	N-His
Predicted MW:	11.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 8 M urea
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001966
Locus ID:	2026
UniProt ID:	P09104 , Q6FHV6
RefSeq Size:	2423
Cytogenetics:	12p13.31
RefSeq ORF:	1302
Synonyms:	HEL-S-279; NSE




[View online »](#)

Summary:

This gene encodes one of the three enolase isoenzymes found in mammals. This isoenzyme, a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates. [provided by RefSeq, Jul 2008]

Protein Pathways:

Glycolysis / Gluconeogenesis, Metabolic pathways, RNA degradation

Product images:

Purified recombinant protein ENO2 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.