

Product datasheet for TP762497

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ATP7B (NM_000053) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human ATPase, Cu++ transporting, beta polypeptide (ATP7B),

transcript variant 1, 200His-482Thr, with N-terminal His tag, expressed in E.coli, 50ug

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

A DNA sequence encoding the region (200His-482Thr) of ATP7B

Tag: N-His

Predicted MW: 31.9kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 10% glycerol

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C after receiving vials.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 000044

Locus ID: 540

UniProt ID: <u>P35670</u>, <u>A0A024RDX3</u>, <u>B7ZLR4</u>

RefSeq Size: 6644

Cytogenetics: 13q14.3

RefSeq ORF: 4395

Synonyms: PWD; WC1; WD; WND






Summary:

This gene is a member of the P-type cation transport ATPase family and encodes a protein with several membrane-spanning domains, an ATPase consensus sequence, a hinge domain, a phosphorylation site, and at least 2 putative copper-binding sites. This protein is a monomer, and functions as a copper-transporting ATPase which exports copper out of the cells, such as the efflux of hepatic copper into the bile. Alternate transcriptional splice variants, encoding different isoforms with distinct cellular localizations, have been characterized. Mutations in this gene have been associated with Wilson disease which is characterized by copper accumulation. [provided by RefSeq, Dec 2019]

Protein Families:

Druggable Genome, Transmembrane

Product images:



Purified recombinant protein ATP7B was analyzed by SDS-PAGE gel and Coomossie Blue Staining.