

## Product datasheet for **TP780003**

### Enterokinase Recombinant Protein (Enterokinase) Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Enterokinase, with C-terminal His tag, expressed in Pichia Pastoris, 500 units
Expression Host:	Other
Tag:	C-His
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS, 50% glycerol
Bioactivity:	1 unit is defined as the amount of enzyme required to cleave 25 µg of a substrate protein to 95% completion in 16 hours at 25°C.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at +4°C upon receipt.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

#### Product images:



[View online »](#)