

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TP790047

Dendra2 Escherichia coli Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified fluorescent protein Dendra2 (photoswitchable), with N-terminal HIS tag, expressed in E. coli, 250ug
Species:	Escherichia coli
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MRGSGSMNTP GINLIKEDMR VKVHMEGNVN GHAFVIEGEG KGKPYEGTQT ANLTVKEGAP LPFSYDILTT AVHYGNRVFT KYPEDIPDYF KQSFPEGYSW ERTMTFEDKG ICTIRSDISL EGDCFFQNVR FKGTNFPPNG PVMQKKTLKW EPSTEKLHVR DGLLVGNINM ALLLEGGGHY LCDFKTTYKA KKVVQLPDAH FVDHRIEILG NDSDYNKVKL YEHAVARYSP LPSQVWSGDS GVYK
Tag:	N-His
Predicted MW:	27.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS, pH 7.4, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>ADE48850</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:

116 — 66 — 45 — 35 — 25 — 18 — 14 —

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US