

## Product datasheet for **TP790049**

### GTGFP Escherichia coli Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified fluorescent protein TGFP, with N-terminal HIS tag, expressed in E. coli, 250ug
Species:	Escherichia coli
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MESDESGLPA MEIECRITGT LNGVEFELVG GGEGTPEQGR MTNKMKSTKG ALTFSPYLLS HVMGYGFYHF GTYPSGYENP FLHAINNGGY TNTRIEKYED GGVLHVSFSY RYEAGRVIGD FKVMGTGFPE DSVIFTDKII RSNATVEHLH PMGDNDLDGS FTRTFSLRDG GYSSVWDSH MHFKSAIHPS ILQNGGPMFA FRRVEEDHSN TELGIVEYQH AFKTPDADAG EE
Tag:	N-His
Predicted MW:	25.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS, pH 7.4, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">ADD23343</a>



[View online »](#)

**Product images:**

