


## Product datasheet for **TP790052**

### CAS9 Streptococcus Pyogenes Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Bacteria <i>S. pyogenes</i> Cas9 endonuclease (Cas9), with N-terminal HIS tag, expressed in <i>E. coli</i> , 50ug
Species:	<i>Streptococcus Pyogenes</i>
Expression Host:	<i>E. coli</i>
Tag:	N-His
Predicted MW:	160 kDa
Concentration:	>0.05 µg/µL as determined by Bradford protein assay
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	10 mM Tris-HCl, pH 7.4, 300 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 50% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -20°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">AE004092</a>

#### Product images:



[View online >](#)