

Product datasheet for **TP790055**

Purified native IgM from human serum, 100 ug Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified native IgM from human serum, 100 ug
Species:	Human
Expression Host:	human serum
Expression cDNA Clone or AA Sequence:	Native form of human IgM
Tag:	Tag Free
Concentration:	>0.5mg/mL as determined by microplate BCA method
Purity:	> 99% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	1 x PBS, pH 7.4, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)