

## Product datasheet for **TP790122**

### p24 HIV Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of HIV-1 p24 Protein, with C-terminal His tag, expressed in E.coli, 50ug
Species:	HIV
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	PIVQNLQGQM VHQAI SPRTL NAWVKVVEEK AFSPEVIPMF SALSEGATPQ DLN TMLNTVG GHQAAMQMLK ETINEEAAEW DRLHPVHAGP IAPGQMREPR GSDIAGTTST LQEQIGWMTH NPPIPVGEIY KRWILGLNK IVRMYSPTSI LDIRQGPKEP FRDYVDRFYK TLRAEQASQE VKNWMTETLL VQNANPDCKT ILKALGPGAT LEEMMTACQG VGGPGHKARV L
Tag:	C-His
Predicted MW:	25.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 90% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS, pH 7.4, 1 mM DTT, 10% glycerol
Endotoxin:	< 1 EU per 1 µg of the protein by the LAL
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">ACI05538</a>



[View online »](#)

**Product images:**

