

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TP790187**

#### BirA Escherichia coli Recombinant Protein


### **Product data:**

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein Biotin ligase (BirA) from Escherichia coli, Tag free, expressed in E.coli
Species:	Escherichia coli
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MKDNTVPLKL IALLANGEFH SGEQLGETLG MSRAAINKHI QTLRDWGVDV FTVPGKGYSL PEPIQLLNAE KILSQLDDGS VAVLPVIDST NQYLLDRIGE LKSGDACVAE YQHAGRGRRG RKWFSPFGAN LYLSMFWRLE QGPAAAIGLS LVIGIVMAEV LRKLGADKVR VKWPNDLYLQ DRKLAGILVE LTGKTGDAAQ IVIGAGINMA MRRVEESVVN QGWITLQEAG INLDRNTLAA MLIRELRAAL ELFEQEGLAP YLSRWEKLDN FINRPVKLII GDKEIFGISR GIDKQGALLL EQDGIIKPWM GGEISLRSAE K
Tag:	Tag Free
Predicted MW:	35.2KDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 6.8, 50 mM NaCl, 5 mM mercaptoethanol, 5% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -20°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

## **Product images:**



Purified recombinant protein BirA was analyzed by SDS-PAGE gel and Coomossie Blue Staining.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US