


Product datasheet for **TP790189**

SARS-CoV-2 Nucleocapsid Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant SARS coronavirus 2 Nucleocapsid Protein (N protein), with C-terminal MYC/DDK tag, expressed in HEK293 cells, 20ug
Species:	SARS-CoV-2
Expression Host:	HEK293
Tag:	Myc-DDK
Predicted MW:	47 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>QHD43423</u>

Product images:



[View online »](#)