

Product datasheet for **TP790195**


LOC113797073 (XM_027347389) Dermatophagoides pteronyssinus Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of troponin C-like, full length, with N-terminal His tag, expressed in E.coli, 50ug
Species:	Dermatophagoides pteronyssinus
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MSVEELTKEQ VQMLRKAFDM FDREKKGHIH TNMVSTILRT LGQTFEEKDL QTLITEIDAD GSGELEFDEF LTLTARFLVE EDAEAMQEEL REAFRMYDKE GNGYIPTSAL REILRALDDK LTNDELDEMI AEIDTDGSGT VDFDEFMEMM TGD
Tag:	N-His
Predicted MW:	19.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 85% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 500 mM NaCl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	XP_027203190



[View online »](#)

Product images:

Purified recombinant protein LOC113797073 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.