

Product datasheet for TP790202


Blot 10 *Blomia tropicalis* Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of group 10 allergen Blo t 10, full length 267Gln, with N-terminal His tag, expressed in E.coli, 50ug
Species:	<i>Blomia tropicalis</i>
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<pre> 1 MGSSHHHHHH SGLVPRGSH MEAIKKKMQA MKLEKDNAID RAEIAEQKSR DANLRAEKSE 61 EEVRALQKKI QQIENELDQV QESLTQANTK LEEKEKSLQT AEGDVAALNR RIQLIEEDLE 121 RSEERLKVAT AKLEEASHSA DESERMKML EHRSITDEER MDGLESQKLE ARMMMAEDADR 181 KYDEVARKLA MVEADLERAE ERAETGETKI VELEEELRVV GNNLKSLEVS EEKAQQREEA 241 YEQQIRMMTG KLKEAEARAE FAERSVQKLQ KEVDRLEDEL VHEKEYKSI SDELDTQFAE 301 LTGY </pre>
Tag:	N-His
Predicted MW:	34.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 90% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 500 mM NaCl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
Locus ID:	ABU97466



[View online »](#)

Product images:

Purified recombinant protein Blot 10 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.