

Product datasheet for **TP790207**


Derf 31 Dermatophagoides farinae Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Der f 31 allergen, SSSLVPRGSHM-17Gly-End, with N-terminal His tag, expressed in E.coli, 50ug
Species:	Dermatophagoides farinae
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	1 MASGVTVATE AKTLYEEVKK DKKYRYIIYH IKDERVIEVE TTGPRDATYS DFVAKLQDYK 61 NECRYCVFDF PANIPVEGGG EKSNSVDRL VLMTWCPESA KIKQKMLYSS SYDALKKALV 121 GVKYVQACD YEEASEEAIA EAFRKGGK
Tag:	N-His
Predicted MW:	16.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 90% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 500 mM NaCl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
Locus ID:	AIO08870



[View online »](#)

Product images:

Purified recombinant protein Derf 31 was analyzed by SDS-PAGE gel and Coomossie Blue Staining.