

Product datasheet for **TR300129**

ZNF568 Human shRNA Plasmid Kit (Locus ID 374900)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ZNF568 Human shRNA Plasmid Kit (Locus ID 374900)
Locus ID:	374900
Synonyms:	ZFP568
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ZNF568 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 374900). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001204835 , NM_001204836 , NM_001204837 , NM_001204838 , NM_001204839 , NM_198539 , NM_198539.1 , NM_198539.2 , NM_198539.3 , NM_001204839.1 , NM_001204836.1 , NM_001204837.1 , NM_001204838.1 , NM_001204835.1 , BC016334 , BC031218 , BC041927 , BC108696 , BC111777 , BC131547 , NM_001204835.2 , NM_001204836.2 , NM_001204837.2
UniProt ID:	Q3ZCX4
Summary:	Has transcriptional repression activity, partially through the recruitment of the corepressor TRIM28 but has also repression activity independently of this interaction. Essential during embryonic development, where it acts as direct repressor of a placental-specific transcript of IGF2 in early development and regulates convergent extension movements required for axis elongation and tissue morphogenesis in all germ layers. Also important for normal morphogenesis of extraembryonic tissues including the yolk sac, extraembryonic mesoderm and placenta. May enhance proliferation or maintenance of neural stem cells. [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).