


## Product datasheet for **TR30013**

### Scrambled shRNA control in pGFP-V-RS shRNA Vector

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Scrambled shRNA control in pGFP-V-RS shRNA Vector
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .
Storage:	Store the vector DNA at -20°C.
Stability:	Vector DNA is stable for at least one year from date of receipt when stored at -20°C.

#### Product images:



User supplied image.

*Alyssa P. for TR30013 on 06/02/2023:*  
E12.5 brain electroplated with scrambled shRNA  
immunostained for Rtkn (purple), Tbr2 (yellow)  
and DAPI (light blue)



[View online »](#)