

## Product datasheet for **TR30016**

### Positive shRNA control against tGFP (pRFP-C-RS Vector)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Positive shRNA control against tGFP (pRFP-C-RS Vector)
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .
Storage:	Store the vector DNA at -20°C.
Stability:	Vector DNA is stable for at least one year from date of receipt when stored at -20°C.



[View online »](#)