

Product datasheet for **TR30016**

Positive shRNA control against tGFP (pRFP-C-RS Vector)

Product data:

Product Type: shRNA Plasmids

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Storage: Store the vector DNA at -20°C.

Stability: Vector DNA is stable for at least one year from date of receipt when stored at -20°C.

