

Product datasheet for **TR30021**

Scrambled shRNA control in pGFP-C-shLenti shRNA Vector

Product data:

Product Type:	shRNA Plasmids
Product Name:	Scrambled shRNA control in pGFP-C-shLenti shRNA Vector
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Storage:	Store the vector DNA at -20°C.
Stability:	Vector DNA is stable for at least one year from date of receipt when stored at -20°C.



[View online »](#)