

Product datasheet for **TR300287**

ZNF189 Human shRNA Plasmid Kit (Locus ID 7743)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ZNF189 Human shRNA Plasmid Kit (Locus ID 7743)
Locus ID:	7743
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ZNF189 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 7743). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<u>NM_001278231</u> , <u>NM_001278232</u> , <u>NM_001278240</u> , <u>NM_003452</u> , <u>NM_197977</u> , <u>NR_103480</u> , <u>NM_003452.1</u> , <u>NM_003452.2</u> , <u>NM_003452.3</u> , <u>NM_197977.1</u> , <u>NM_197977.2</u> , <u>NM_001278240.1</u> , <u>NM_001278232.1</u> , <u>NM_001278231.1</u> , <u>BC092425</u> , <u>BC065718</u> , <u>BC126120</u> , <u>BC143876</u> , <u>BC143877</u> , <u>BM820380</u> , <u>BM996467</u> , <u>NM_003452.4</u> , <u>NM_001278231.2</u> , <u>NM_001278232.2</u>
UniProt ID:	<u>O75820</u>
Summary:	Kruppel-like zinc finger proteins such as ZNF189 contain a conserved stretch of 7 amino acids that connects a variable number of DNA-binding zinc finger repeats of the cys(2)his(2) (C2H2) type (summarized by Odeberg et al., 1998 [PubMed 9653648]). Approximately 30% of human Kruppel-like zinc finger proteins contain an N-terminal Kruppel-associated box (KRAB) domain. The KRAB domain consists of approximately 75 amino acids that may be subdivided into an A box, which is present in every KRAB domain and is essential for transcriptional repression, and a B box, which is not always present.[supplied by OMIM, May 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).