

Product datasheet for **TR30031**

Scrambled shRNA control in pRFP-C-shLenti vector

Product data:

Product Type:	shRNA Vectors
Control Type:	Negative
Fluorescent Label:	None
E. coli Selection:	Chloramphenicol
Mammalian Cell Selection:	Puromycin
Viral Packaging:	Lentiviral
Feature:	Transfection grade purified shRNA plasmid carrying a 29-mer scrambled sequence cassette. This product is a negative control for OriGene's gene specific HuSH shRNA lentiviral products.
Vector Sequence:	Download
Storage:	Store the vector DNA at -20°C.
Stability:	Vector DNA is stable for at least one year from date of receipt when stored at -20°C.



[View online »](#)